


# Codex<sup>®</sup>

## HHDH Panel

Screening Protocol

Document # PRO-015-03

### Reactions of Interest



### Codex<sup>®</sup> HHDH Panel Plate Information

- Recommended storage temperature for the plates is -20 °C.
- The plates are from Corning (product number 3961) and can be purchased from VWR (catalog number 29445-166).

### Codex<sup>®</sup> HHDH Panel Materials and Reagents

- Two 96-well plates, provided.
- DMSO, not provided.
- 100 mM Tris-SO<sub>4</sub> buffer, pH 7, not provided.

## Codex® HHDH Panel

### Codex® HHDH Panel Screening Procedure

1. Defrost plates, which contain 100  $\mu$ L lysate per well, at 4 °C. Centrifuge (4000 rpm, 2 min, 4 °C).
2. Assay setup (final reaction volume 300  $\mu$ L per well):
  - a. Unseal the plates.
  - b. Add the substrate solution to each well.
  - c. Seal the assay plate (180 °C, 3 sec).
3. Shake the plates on titer plate shaker (low speed) at room temperature for 24 hours.

| Reagent Concentration for Epoxide formation                          | Volume per well | Volume per plate | Final concentration                              |
|----------------------------------------------------------------------|-----------------|------------------|--------------------------------------------------|
| <b>Lysate</b>                                                        | 100 $\mu$ L     | 9.60 mL          | 40% v/v                                          |
| <b>Substrate solution</b><br>75 mM halohydrin<br>8.3% DMSO<br>Buffer | 150 $\mu$ L     | 14.4 mL          | 60% v/v<br>45 mM halohydrin<br>5% DMSO<br>Buffer |
| <b>Total Volume</b>                                                  | 250 $\mu$ L     | 24.0 mL          |                                                  |

| Reagent Concentration for Epoxide formation                                                  | Volume per well | Volume per plate | Final concentration                                                      |
|----------------------------------------------------------------------------------------------|-----------------|------------------|--------------------------------------------------------------------------|
| <b>Lysate</b>                                                                                | 100 $\mu$ L     | 9.60 mL          | 40% v/v                                                                  |
| <b>Substrate solution</b><br>10 mM epoxide<br>20 mM nucleophile<br>Up to 8.3% DMSO<br>Buffer | 150 $\mu$ L     | 14.4 mL          | 60% v/v<br>10 mM epoxide<br>20 mM nucleophile<br>Up to 5% DMSO<br>Buffer |
| <b>Total Volume</b>                                                                          | 250 $\mu$ L     | 24.0 mL          |                                                                          |

## Codex® HHDH Panel

### Codex® HHDH Panel Work-Up and Analysis

1. Reaction quenching:
  - a. Centrifuge the assay plates (4000 rpm, 1 min, 4 °C).
  - b. Unseal the assay plates.
  - c. Add 1000 µL of quenching solvent (e.g. acetonitrile, methanol or MTBE) to each well.
  - d. Seal the quenched plates (180 °C, 3 sec).
2. Shake the plates on titer plate shaker (medium-high speed) at room temperature for 20 minutes.
3. HPLC sample preparation:
  - a. Centrifuge the assay plates (4000 rpm, 20 min, 4 °C).
  - b. Unseal the quenched plates.
  - c. Transfer 200 µL of samples from quenched plates to round-bottom shallow well plates (Corning 3365).
4. Seal the shallow well plates (180 °C, 2 sec). Analyze using HPLC method of choice.

For further information or any questions,  
please contact us at [sales@codexis.com](mailto:sales@codexis.com)